


TDS 接线:



④ Apply

- 清零:
- ① To monitor
 - ② Proof limit
 - ③ Initial-in + Start

激活通道 MEAs

- ① Simple measurement setting
- ② not use simple \Rightarrow Enter
- ③ To monitor
- ④ Monitor channel setting

④ Monitor channel setting

PC端设置

Measurement \rightarrow View measurement panel (第一次要保存)
 \rightarrow 下一步 \rightarrow 取消 Simple measure \rightarrow \checkmark Initial 0-9 \rightarrow normal
Set monitor \rightarrow \checkmark 选全部 \rightarrow Insert

时间道路设置

CH10 = CH(001) 单位, 小数点
时间设定要比当下时间提前